



NANO

1.155

 **Allowed user weight**
130 kg

 **Suitable for passenger transport**
Tested according to ISO 7176-19 (max. user weight 130 kg)





Vario element



Shapely frame concept



Steering wheel connection



Standard with folding back

Pure design

EFFICIENCY

- Permanent provision with adjustment possibilities to functional changes reduces follow-up costs
- Extensive serial equipment minimises necessity for expensive customised parts
- Low wear due to high quality attachments

THERAPEUTIC BENEFIT

- Perfect easy running and little own weight support the independence of the user
- Adaptable to functional changes and individual needs of the user through adjustment of the back, seat and chassis
- Extensive equipment assortment for long-term users

FIELD OF USE

- For permanent and long term provision
- For paraplegic provision
- Very good suitability for many disease patterns, even for severely encumbered
- Users without function of the lower extremities (lateral transfer, driven through the handrims)

ORDER

COST ESTIMATE

CUSTOM DESIGN

Please fax the completed order form to: 05733 922 9311

Customer service: 05733 922 311

Billing address / Customer ref. no.:

Delivery address:

Consignment / Comments:

Order quantity:

MEYRA



Allowed user weight

130 kg



Suitable for passenger transport

Tested according to ISO 7176-19 (max. user weight 130 kg)



Hints

Our Terms and Conditions as stated at www.meyra.de/AGB shall apply
You can find information on the measurement process at www.meyra.de/messverfahren

Valid from 01/01/2021 - 12/31/2021.
Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check



Technical specifications

Seat widths in mm	300/320/340/360/380/400/420/440/460/480/500
Ready for use width in mm	Seat width + 170
Seat depth in mm	380/400/420/440/460/480/500
Seat height front/rear in mm	450 to 530/360 to 510
Overall height in mm	Rear seat height + back support height
Back support height in mm	variable from 250 to 450 +/- 15
Length with footrest in mm	720 - 920 (depending on frame length, wheel size and axle position)
User weight	130 kg
Empty weight	9,5 kg
Permissible overall weight	145 kg
Transport weight	6,5 kg
Driving wheel	600 mm(24 inch) / 635 mm(25 inch) / 660 mm(26 inch)
Brake test according to STANDARD	7°

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

00 Seat width / Seat depth

Seat width/Seat depth	Short frame			Medium frame		Long frame	
	SD 380 mm	SD 400 mm	SD 420 mm	SD 440 mm	SD 460 mm	SD 480 mm	SD 500 mm
SW 300 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 320 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 340 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 360 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 380 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 400 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 420 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 440 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 460 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 480 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 500 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

including utensil pocket under the seat belt

01 Frame forms

CODE	DESCRIPTION	ANGLE
<input type="checkbox"/> 346	Standard front frame	80°
<input type="checkbox"/> 347	Compact front frame ⚠ Only in combination with frame reduction CODE 4950	88°

02 Frame heeling

CODE	DESCRIPTION
<input type="checkbox"/> 4950	With frame heeling (20 mm each side)
<input type="checkbox"/> 4959	Without frame heeling

Please select by ticking the relevant boxes

S *SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

03 Castor wheel (without seat cushion) in combination with front seat height (FSH)

Code/Art.	Castor	Front seat height in mm (FSH)								
		450	460	470	480	490	500	510	520	530
151	Ø 100 Solid rubber	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
325	Ø 125 Solid rubber		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
327/908	Ø 142 Soft-solid rubber black			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
327/909	Ø 142 Soft-solid rubber grey			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

X-Trend castor wheel (without seat cushion) in combination with front seat height (FSH)

Code/Art.	Castor	Front seat height in mm (FSH)								
		450	460	470	480	490	500	510	520	530
1068963	Ø 100 X-Trend Lightwheel 1068963	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1072385	Ø 100 X-Trend Lightwheel 1072385	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1072386	Ø 125 X-Trend Lightwheel 1072386		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Aluminum castor wheel (without seat cushion) in combination with front seat height (FSH)

FSH (mm)	450	460	470	480	490	500	510	520
A	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
B		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Code	Castor	Rim color					
		Silver	black	Green	Red	Blue	Gold
689	Ø 100 x 37 mm, Aluminum rim, Solid rubber	<input type="checkbox"/> A	<input type="checkbox"/> A	<input type="checkbox"/> A	<input type="checkbox"/> A	<input type="checkbox"/> A	<input type="checkbox"/> A
689	Ø 125 x 37 mm, Aluminum rim, Solid rubber	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B

04 Seat cushions

CODE	DESCRIPTION
<input type="checkbox"/> 560	Seat cushion 30 mm
<input type="checkbox"/> 561	Seat cushion 60 mm

Please select by ticking the relevant boxes

S *SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

05 Drive wheels in combination with seat height

RSH (mm)	360	370	380	390	400	410	420	430	440	450	460	470	480	490	500	510
A	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
B			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
C					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

CODE xR	RSH	DESCRIPTION	Ø
<input type="checkbox"/> 651	A	Standard wheel, incl. quick release axle (pair)	24"
<input type="checkbox"/> 2651	A	Standard wheel black, incl. quick release axle (pair)	24"
<input type="checkbox"/> 848	A	SILVERDESIGN with radial spokes and high flange rim incl. quick-release axle (pair)	24"
<input type="checkbox"/> 752	A	Lightwheel incl. RightRun Light tyres and aluminium handrim black ! Only in combination with 0° wheel camber CODE 423. Recommended up to 90 kg user weight	24"
<input type="checkbox"/> 695	A	Active-wheel with radial spokes, high flange rim, hollow rim (Jumbo rim) in black ! Recommended from 100 kg user weight or in case of dynamic driving habits	24"
<input type="checkbox"/> 695	B	Active-wheel with radial spokes, high flange rim, hollow rim (Jumbo rim) in black ! Recommended from 100 kg user weight or in case of dynamic driving habits	25"
<input type="checkbox"/> 695	C	Active-wheel with radial spokes, high flange rim, hollow rim (Jumbo rim) in black ! Recommended from 100 kg user weight or in case of dynamic driving habits	26"

Drive wheels X-Trends

RSH (mm)	360	370	380	390	400	410	420	430	440	450	460	470	480
A	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Item	RSH	DESCRIPTION	Ø
<input type="checkbox"/> 1081354	A	Active wheel with radial spoke, Schwalbe ONE tires, hub, rim and handrims BLUE	24x1"
<input type="checkbox"/> 1073640	A	Day-Lite drive wheel black, including Speedlite tires and aluminum handrims black	25x1"
<input type="checkbox"/> 1081356	A	Active wheel with radial spoke, Schwalbe ONE tires, hub, rim and handrims RED	24x1"
<input type="checkbox"/> 1081355	A	Active wheel with radial spoke, Schwalbe ONE tires, hub, rim and handrims GOLD	24x1"
<input type="checkbox"/> 1068713	A	Black Light incl. Marathon-Plus-Evolution and aluminium handrim	24x1"
<input type="checkbox"/> 1073640	A	Spinergy-wheel 12 spokes black, incl. Marathon-Plus-Evolution and stainless steel handrim ! White or orange spokes possible on request!	24x1"

Horizontal drive wheel position (balance point setting)


Position	PASSIV (Pos./mm)			MEDIUM (Pos./mm)				ACTIVE (Pos./mm)**			
	1/50	2/60	3/70	4/80	5/90	6/100	7/110	8/120	9/130	10/140	11/150
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

! Rear wheel position 8, 9, 10 und 11 not in combination with frame short CODE 351

Please select by ticking the relevant boxes *SBI - Sonderbau/Individuallösung Standard version: If no other selection is made, the standard version will be delivered.**06 Drive wheel options**

CODE	DESCRIPTION
<input type="checkbox"/> 880	Wheelbase extension 0 or 735 mm

Wheel camber







CODE	DESCRIPTION
<input type="checkbox"/> 423	0°-wheel camber (track adjustable)
<input type="checkbox"/> 447	3°-wheel camber (track adjustable)
<input type="checkbox"/> 454	6°-wheel camber (track adjustable)  not in combination with light brake CODE 2595 max. user weight 100 kg

Please select by ticking the relevant boxes


*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

07 Hand rims (not in combination with Xtrends wheels)

CODE	DESCRIPTION	Ø
<input type="checkbox"/> 493	Stainless steel, polished	24"/25"/26"
<input type="checkbox"/> 494	Aluminium  Not in combination with Active-Wheel CODE 695	24"
<input type="checkbox"/> 445	Titanium, with 6 brackets  Only in combination with active wheel CODE 695 or Spinergy wheel ART 1073640	24"/25"/26"
<input type="checkbox"/> 168	Optigripp-special handrim for tetraplegics	24"
Curve Hand rim in conjunction with Active Wheel CODE 695 or standard wheel CODE 651 and 2651		
<input type="checkbox"/> 849	Curve L, aluminium black anodised, ergonomic handrim profile, H 31 mm  only in conjunction with Active Wheel CODE 695 or standard wheel CODE 651 and 2651	24"
<input type="checkbox"/> 850	Curve L Tetra Grip, anti-slip coated, aluminum, suitable for quadriplegics, ergonomic handrim profile, H 31 mm  only in conjunction with Active Wheel CODE 695 or standard wheel CODE 651 and 2651	24"
<input type="checkbox"/> 853	Curve GEKKO, aluminum black anodized with exchangeable rubber lip, ergonomic handrim profile, H 31 mm  only in conjunction with Active Wheel CODE 695 or standard wheel CODE 651 and 2651	24"
Curve Hand rim in conjunction with Spinergy wheel article number: 1073640		
<input type="checkbox"/> 1078827	Curve L, aluminium black anodised, ergonomic handrim profile, H 31 mm  only in conjunction with Spinergy wheel article number: 1073640	24"
<input type="checkbox"/> 1078828	Curve L Tetra Grip, anti-slip coated, aluminum, suitable for quadriplegics, ergonomic handrim profile, H 31 mm  only in conjunction with Spinergy wheel article number: 1073640	24"
<input type="checkbox"/> 1080439	Curve GEKKO, aluminum black anodized with exchangeable rubber lip, ergonomic handrim profile, H 31 mm  only in conjunction with Spinergy wheel article number: 1073640	24"

Hand rim options

CODE	DESCRIPTION	
<input type="checkbox"/> 166	Handrim coating made of silicone, smooth, black  Not in combination with Active-Wheel CODE 695	24"
<input type="checkbox"/> 692	Handrim distance narrow	24"
<input type="checkbox"/> 693	Handrim distance far	24"

Please select by ticking the relevant boxes


*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

08 Tyres (not combinable with X-Trends wheels)

CODE	DESCRIPTION	Ø		
<input type="checkbox"/> 496	Pneumatic tires 8 bar	24x1"		
<input type="checkbox"/> 897	Schwalbe Marathon Plus Evolution tires, puncture resistant	24x1"	25x1"	26x1"
<input type="checkbox"/> 363	Schwalbe ONE tires, no profile, puncture resistant	24x1"		
<input type="checkbox"/> 395	Maxxis VN Master Alpha Protect tire black, puncture resistant, 10 bar	24x1"		

09 Brakes

CODE	DESCRIPTION
<input type="checkbox"/> 305	Pressure brake for user
<input type="checkbox"/> 2595	Light pressure brake made of aluminum, swing away
<input type="checkbox"/> 681	Brake lever extension, foldable  only in combination with CODE 305 pressure brake

10 Back strap height (BH)

DESCRIPTION	Back strap height in mm				
	250	300	350	400	450
Back strap, adaptable via Velcro straps adjustable +/- 15 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

11 Back support

CODE	DESCRIPTION
<input type="checkbox"/> 410	Back support with push handle short, 50 mm height adjustable
<input type="checkbox"/> 411	Back support without push handle
<input type="checkbox"/> 502	Back support with push handle mounted to the frame tube, continuously height adjustable to 150 mm
<input type="checkbox"/> 141	Push handles, swivel-away

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

11 Backrest options

CODE	DESCRIPTION
<input type="checkbox"/> 2032	ErgoSeat Pro thorax pelottes, small, pair ⚠ Recommended at back height 300 - 400 mm
<input type="checkbox"/> 2021	ErgoSeat Pro thorax pelottes, large, pair ⚠ Recommended at back height 400 - 450 mm

12 Arm supports and arm pads

CODE	DESCRIPTION	Functions
<input type="checkbox"/> 748	Aluminum side guard in frame color with cold protection insert including height-adjustable armrest (screwed)	1,5,6
<input type="checkbox"/> 100	Side guard made of aluminum in frame color with cold protection insert, screwed tight ⚠ i.C. with wheel size 24 x 1 inch from RSH 440mm i.C. with wheel size 25 x 1 inch from RSH 470mm not i.C. with wheel size 26 x 1 inch	3,4,6
<input type="checkbox"/> 624/100	Side guard made of aluminum in frame color with cold protection insert, screwed tight ⚠ i.C. with wheel size 24 x 1 inch from RSH 390mm-430mm i.C. with wheel size 25 x 1 inch from RSH 420mm-460mm i.C. with wheel size 26 x 1 inch from RSH 450mm-510mm	3,4,6
<input type="checkbox"/> 2624/100	Side guard made of aluminum in frame color with cold protection insert, screwed tight ⚠ i.C. with wheel size 24 x 1 inch to HSH 380mm i.C. with wheel size 25 x 1 inch to HSH 410mm i.C. with wheel size 26 x 1 inch to HSH 440mm	3,4,6
<input type="checkbox"/> 783	Carbon arm support	
<input type="checkbox"/> 11	Side guard removable ⚠ only in combination with CODE 100, 783 oder 748	
<input type="checkbox"/> 24	Swing-away side panel bracket ⚠ only in combination with CODE 100, 783 or without side panel.	

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

13 Leg supports and footrests

LSL (mm) **360** 370 380 390 400 410 420 430 440 450 460 470 480

CODE xR	LSL (mm)	DESCRIPTION
<input type="checkbox"/> 76	<input checked="" type="checkbox"/> A	Swiveling footrest
<input type="checkbox"/> 311	<input checked="" type="checkbox"/> A	Tubular footrest, aluminium, height-adjustable, black matt
<input type="checkbox"/> 312	<input checked="" type="checkbox"/> A	Footrest one-piece, height / depth and angle adjustable
<input type="checkbox"/> 4312	<input checked="" type="checkbox"/> A	LED footrest one-piece, made of strong aluminum in frame color, height / depth and angle adjustable, incl. Bluetooth module ⚠ i.C. with CODE 4959 from SW 360mm i.C. with CODE 4950 from SW 400mm
<input type="checkbox"/> 4313	<input checked="" type="checkbox"/> A	LED footrest one-piece, made of strong aluminum in frame color, height / depth and angle adjustable ⚠ i.C. with CODE 4959 from SW 360mm i.C. with CODE 4950 from SW 400mm
<input type="checkbox"/> 359	Calf belt	

Please select by ticking the relevant boxes

S *SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

14 Frame colours

CODE	DESCRIPTION
<input type="checkbox"/> 130	Titanium grey
<input type="checkbox"/> 183	Metallic silver
<input type="checkbox"/> 186	Magic Blue
<input type="checkbox"/> 216	Deep black
<input type="checkbox"/> 206	Matt black
<input type="checkbox"/> 224	Blue-Metallic
<input type="checkbox"/> 227	Red Metallic
<input type="checkbox"/> 4023	Signal white
<input type="checkbox"/> 4040	Downtown Gray
<input type="checkbox"/> 4041	Magicblack Matt
<input type="checkbox"/> 4042	Blue Metallic Matt
<input type="checkbox"/> 4043	Brown metallic matt
<input type="checkbox"/> 4044	Red Metallic Matt
<input type="checkbox"/> 4045	Signal white matt
<input type="checkbox"/> 4046	Green Metallic Matt
<input type="checkbox"/> 4047	Purple metallic matt
<input type="checkbox"/> 4048	Copper metallic
<input type="checkbox"/> 4913	Cross orange
<input type="checkbox"/> 4922	Metallic green

Other frame colors are available on request for an extra charge!

Color of choice

Colour code **130 183 186 206 216 224 227 4023 4040 4041 4042 4043 4044 4045 4046 4047 4048 4913 4922**

F

CODE xR	Colour code	DESCRIPTION
<input type="checkbox"/> 2205	F	side guards/footrest in custom colour ! only in combination with Code 100 side guard fixed

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

15 Back strap colours




CODE	DESCRIPTION
<input type="checkbox"/> 235	Trevira black
<input type="checkbox"/> 292	Net fabric, breathable , black

16 accessories

Hand and spoke guard

CODE	DESCRIPTION
<input type="checkbox"/> 489	Transparent

Support castors


CODE	DESCRIPTION
<input type="checkbox"/> 728	Support castors, 1 pair, swivelling (recommended)  not in combination with step tube CODE 309
<input type="checkbox"/> 729	1 support castor, swivelling, mounted left  not in combination with the display tube CODE 309/4930
<input type="checkbox"/> 730	1 support castor, swivelling, mounted right  not in combination with entry pipe CODE 309/4931

Pedal tubes

CODE	DESCRIPTION
<input type="checkbox"/> 309	Step tubes (1 pair)  not in combination with support rollers CODE 728, CODE 729 or CODE 730
<input type="checkbox"/> 4930	1 step tube, mounted on the left  not in combination with support rollers CODE 728 or CODE 729
<input type="checkbox"/> 4931	1 step tube, mounted on the right  not in combination with support rollers CODE 728 or CODE 730

Please select by ticking the relevant boxes *SBI - Sonderbau/Individuallösung Standard version: If no other selection is made, the standard version will be delivered.

Accessories

CODE	DESCRIPTION
<input type="checkbox"/> 833	Retaining strap with buckle
<input type="checkbox"/> 970	Walking aid holder  only in combination with CODE 309/4390/4391
<input type="checkbox"/> 819	Reflectors
<input type="checkbox"/> 951	Accessory net

Remarks

Date

City

Signature