


NANO X

1.156

 **Allowed user weight**
130 kg






Vario element



Shapely frame concept



Castor wheel connection milled from solid material



Double cross brace for rigid frame-like ride performance

Design your Freedom

EFFICIENCY

- Permanent provision with adjustment possibilities to functional changes reduces follow-up costs
- Extensive serial equipment minimises necessity for expensive customised parts
- Low wear due to high quality attachments

THERAPEUTIC BENEFIT

- Perfect easy running and little own weight support the independence of the user
- Adaptable to functional changes and individual needs of the user through adjustment of the back, seat and chassis
- Extensive equipment assortment for long-term users

FIELD OF USE

- For permanent and long term provision
- For paraplegic provision
- Very good suitability for many disease patterns, even for severely encumbered
- Users without function of the lower extremities (lateral transfer, driven through the handrims)

ORDER

COST ESTIMATE

CUSTOM DESIGN

Please fax the completed order form to: 05733 922 9311

Customer service: 05733 922 311

Billing address / Customer ref. no.:

Delivery address:

Consignment / Comments:

Order quantity:

MEYRA



Allowed user weight

130 kg



Hints

Our Terms and Conditions as stated at www.meyra.de/AGB shall apply
You can find information on the measurement process at www.meyra.de/messverfahren

Valid from 01/01/2021 - 12/31/2021.
Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check



Technical specifications

Seat widths in mm	320/340/360/380/400/420/440/460/480
Ready for use width in mm	480/500/520/540/560/580/600/620/640/660/680
Seat depth in mm	380/400/420/440/460/480/500
Seat height front/rear in mm	450 to 530/360 to 510
Overall height in mm	Rear seat height + back support height
Back support height in mm	variable from 250 to 450 +/- 15
Length with footrest in mm	720 - 920 (depending on frame length, wheel size and axle position)
User weight	130 kg
Empty weight	ab 10 kg
Permissible overall weight	145 kg
Transport weight	ab 7 kg
Driving wheel	600 mm(24 inch) / 635 mm(25 inch) / 660 mm(26 inch)
Brake test according to STANDARD	7°

Please select by ticking the relevant boxes

S *SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

00 Seat width / Seat depth

Seat width/Seat depth	Short frame			Medium frame		Long frame	
	SD 380 mm	SD 400 mm	SD 420 mm	SD 440 mm	SD 460 mm	SD 480 mm	SD 500 mm
SW 320 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 340 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 360 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 380 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 400 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 420 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 440 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 460 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 480 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

including utensil pocket under the seat belt

01 Frame forms

CODE	DESCRIPTION	ANGLE
<input type="checkbox"/> 346	Standard front frame	80°
<input type="checkbox"/> 347	Compact front frame ! Only in combination with frame reduction CODE 4950 not i.C. with CODE 327 Ø 142 Soft-solid rubber	88°
<input type="checkbox"/> 721	Force node tie-down system for passenger transport in vehicles	

02 Frame heeling

CODE	DESCRIPTION
<input type="checkbox"/> 4950	With frame heeling (20 mm each side)
<input type="checkbox"/> 4959	Without frame heeling

Please select by ticking the relevant boxes

S *SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

03 Castor wheel (without seat cushion) in combination with front seat height (FSH)

Please select your desired seat height by selecting a possible combination of seat height and wheel size.
The max permissible seat height difference between front and rear seat height is 100 mm!

Recommendation: Solid rubber steering wheels for heavier user weights

Code/Art.	Castor	Front seat height in mm (FSH)								
		450	460	470	480	490	500	510	520	530
151	Ø 100 Solid rubber	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
325	Ø 125 Solid rubber		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
327/908	Ø 142 Soft-solid rubber black		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
327/909	Ø 142 Soft-solid rubber grey			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

X-Trend castor wheel (without seat cushion) in combination with front seat height (FSH)

Code/Art.	Castor	Front seat height in mm (FSH)								
		450	460	470	480	490	500	510	520	530
1068963	Ø 100 X-Trend Lightwheel 1068963	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
1072385	Ø 100 X-Trend Lightwheel 1072385	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
1072386	Ø 125 X-Trend Lightwheel 1072386		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Aluminum castor wheel (without seat cushion) in combination with front seat height (FSH)

FSH (mm)	450	460	470	480	490	500	510	520
A	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
B		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Code	Castor	Rim color					
		Silver	black	Orange	Red	Blue	White
689	Ø 100 x 37 mm, Aluminum rim, Solid rubber	<input type="checkbox"/> A	<input type="checkbox"/> A	<input type="checkbox"/> A	<input type="checkbox"/> A	<input type="checkbox"/> A	<input type="checkbox"/> A
689	Ø 125 x 37 mm, Aluminum rim, Solid rubber	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B

Please select by ticking the relevant boxes

S *SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

04 Seat cushions

CODE	DESCRIPTION
<input type="checkbox"/> 560	Seat cushion 30 mm
<input type="checkbox"/> 561	Seat cushion 60 mm

05 Drive wheels in combination with seat height

*Bitte wählen Sie Ihre gewünschte Sitzhöhe durch Auswahl einer möglichen Kombination aus Sitzhöhe und Radgröße. 25" und 26" siehe Xtrends Antriebsräder

RSH (mm)	360	370	380	390	400	410	420	430	440	450	460	470	480	490	500	510
A	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
B			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
C					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

CODE xR	RSH	DESCRIPTION	Ø
<input type="checkbox"/> 651	A	Standard wheel, incl. quick release axle (pair)	24"
<input type="checkbox"/> 2651	A	Standard wheel black, incl. quick release axle (pair)	24"
<input type="checkbox"/> 706	A	Drum braked wheel with radial spokes and drum brake for accompanying person ⚠ Only in combination with push handles CODE 502 only i.C. with CODE 423 0° wheel camber	24"
<input type="checkbox"/> 848	A	SILVERDESIGN with radial spokes and high flange rim incl. quick-release axle (pair)	24"
<input type="checkbox"/> 752	A	Lightwheel incl. RightRun Light tyres and aluminium handrim black ⚠ Only in combination with 0° wheel camber CODE 423. Recommended up to 90 kg user weight	24"
<input type="checkbox"/> 695	A	Active-wheel with radial spokes, high flange rim, hollow rim (Jumbo rim) in black ⚠ Recommended from 100 kg user weight or in case of dynamic driving habits	24"
<input type="checkbox"/> 695	B	Active-wheel with radial spokes, high flange rim, hollow rim (Jumbo rim) in black ⚠ Recommended from 100 kg user weight or in case of dynamic driving habits	25"
<input type="checkbox"/> 695	C	Active-wheel with radial spokes, high flange rim, hollow rim (Jumbo rim) in black ⚠ Recommended from 100 kg user weight or in case of dynamic driving habits	26"

⚠ The max permissible seat height difference between front and rear seat height is 100 mm!

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

Drive wheels X-Trends

RSH (mm)	360	370	380	390	400	410	420	430	440	450	460	470	480	490	500
A	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
B			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Item	RSH	DESCRIPTION	Ø
<input type="checkbox"/> 1085185	A	Spinergy Carbon Antriebsrad, Speichen schwarz inkl. Schwalbe One Bereifung und Edelstahl Greifreifen	24x1"
<input type="checkbox"/> 1081354	A	Active wheel with radial spoke, Schwalbe ONE tires, hub, rim and handrims BLUE	24x1"
<input type="checkbox"/> 1081356	A	Active wheel with radial spoke, Schwalbe ONE tires, hub, rim and handrims RED	24x1"
<input type="checkbox"/> 1081355	A	Active wheel with radial spoke, Schwalbe ONE tires, hub, rim and handrims GOLD	24x1"
<input type="checkbox"/> 1068713	A	Black Light incl. Marathon-Plus-Evolution and aluminium handrim ⚠ Nicht i.V. mit CODE 701 Kraftknotensatz für den Personentransport im Auto	24x1"
<input type="checkbox"/> 1073640	A	Spinergy-wheel 12 spokes black, incl. Marathon-Plus-Evolution and stainless steel handrim. Upon inquiry white or orange spokes possible	24x1"
<input type="checkbox"/> 1087341	B	Day-Lite Antriebsrad	25x1"

Horizontal drive wheel position (balance point setting)

Messverfahren: Hinterkante Rückenrohr bis Achsmittle in mm

	PASSIV (Pos./mm)		MEDIUM (Pos./mm)		ACTIVE (Pos./mm)**		
	1/45	2/57	3/69	4/81	5/93	6/105	7/145
Position	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

⚠ Antriebsradposition 7 nicht in Verbindung mit Rahmen kurz CODE 351

06 Wheel camber



CODE	DESCRIPTION
<input type="checkbox"/> 423	0° wheel camber
<input type="checkbox"/> 447	3°-wheel camber (track adjustable)

Please select by ticking the relevant boxes


*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

07 Hand rims (not in combination with Xtrends wheels)

CODE	DESCRIPTION	Ø
<input type="checkbox"/> 493	Stainless steel, polished	24"/25"/26"
<input type="checkbox"/> 494	Aluminium  Not in combination with Active-Wheel CODE 695	24"
<input type="checkbox"/> 445	Titanium, with 6 brackets  Only in combination with active wheel CODE 695 or Spinergy wheel ART 1073640	24"/25"/26"
<input type="checkbox"/> 168	Optigripp-special handrim for tetraplegics	24"
Curve Hand rim in conjunction with Active Wheel CODE 695 or standard wheel CODE 651 and 2651		
<input type="checkbox"/> 849	Curve L, aluminium black anodised, ergonomic handrim profile, H 31 mm  only in conjunction with Active Wheel CODE 695 or standard wheel CODE 651 and 2651	24"
<input type="checkbox"/> 850	Curve L Tetra Grip, anti-slip coated, aluminum, suitable for quadriplegics, ergonomic handrim profile, H 31 mm  only in conjunction with Active Wheel CODE 695 or standard wheel CODE 651 and 2651	24"
<input type="checkbox"/> 853	Curve GEKKO, aluminum black anodized with exchangeable rubber lip, ergonomic handrim profile, H 31 mm  only in conjunction with Active Wheel CODE 695 or standard wheel CODE 651 and 2651	24"
Curve Hand rim in conjunction with Spinergy wheel article number: 1073640		
<input type="checkbox"/> 1078827	Curve L, aluminium black anodised, ergonomic handrim profile, H 31 mm  only in conjunction with Spinergy wheel article number: 1073640	24"
<input type="checkbox"/> 1078828	Curve L Tetra Grip, anti-slip coated, aluminum, suitable for quadriplegics, ergonomic handrim profile, H 31 mm  only in conjunction with Spinergy wheel article number: 1073640	24"
<input type="checkbox"/> 1080439	Curve GEKKO, aluminum black anodized with exchangeable rubber lip, ergonomic handrim profile, H 31 mm  only in conjunction with Spinergy wheel article number: 1073640	24"

Hand rim options

CODE	DESCRIPTION	
<input type="checkbox"/> 166	Handrim coating made of silicone, smooth, black  Not in combination with Active-Wheel CODE 695	24"
<input type="checkbox"/> 692	Handrim distance narrow	24"
<input type="checkbox"/> 693	Handrim distance far	24"

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

08 Tyres (not combinable with X-Trends wheels)

CODE	DESCRIPTION	Ø		
<input type="checkbox"/> 496	Pneumatic tires 8 bar	24x1"		
<input type="checkbox"/> 395	Maxxis VN Master Alpha Protect tire black, puncture resistant, 10 bar	24x1"		
<input type="checkbox"/> 897	Schwalbe Marathon Plus Evolution tires, puncture resistant	24x1"	25x1"	26x1"
<input type="checkbox"/> 363	Schwalbe ONE tires, no profile, puncture resistant	24x1"		

⚠ Must be filled if Standard, Silverdesign or Active wheel is chosen

Note: Tyres consist of a rubber compound which may leave permanent marks or be hard to remove from some surfaces (e.g. plastic, wooden or parquet flooring, rugs and carpets).

We cannot accept liability for any abrasion damage to flooring.

09 Brakes

CODE	DESCRIPTION
<input type="checkbox"/> 305	Pressure brake for user
<input type="checkbox"/> 2595	Light pressure brake made of aluminum, swing away
<input type="checkbox"/> 681	Brake lever extension, foldable

⚠ only in combination with CODE 305 pressure brake

10 Back strap height (BH)

DESCRIPTION	Back strap height in mm				
	250	300	350	400	450
Back strap, adaptable via Velcro straps adjustable +/- 15 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

11 Back support



*Standard 90° zum Sitz

Verstellbar im Winkel und verriegelbar

Rückenlehnenwinkel: -5°/ +5°/ + 10°

CODE	DESCRIPTION
<input type="checkbox"/> 410	Back support with push handle short, 50 mm height adjustable
<input type="checkbox"/> 411	Back support without push handle
<input type="checkbox"/> 502	Back support with push handle mounted to the frame tube, continuously height adjustable to 150 mm
<input type="checkbox"/> 141	Push handles, swivel-away

11 Backrest options

CODE	DESCRIPTION
<input type="checkbox"/> 2032	ErgoSeat Pro thorax pelottes, small, pair  Recommended at back height 300 - 400 mm
<input type="checkbox"/> 2021	ErgoSeat Pro thorax pelottes, large, pair  Recommended at back height 400 - 450 mm

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

12 Arm supports and arm pads

Functional properties:

- 1) fixed
- 2) swivel up
- 3) detachable
- 4) lockable
- 5) height-adjustable armrest
- 6) with wheel cover
- 7) space-saving swivel in

CODE	DESCRIPTION	Functions
<input type="checkbox"/> 100	Side guard made of aluminum in frame color with cold protection insert, screwed tight ⚠ i.C. with wheel size 24 x 1 inch from RSH 440mm i.C. with wheel size 25 x 1 inch from RSH 470mm not i.C. with wheel size 26 x 1 inch	1,6
<input type="checkbox"/> 100/624	Side guard made of aluminum in frame color with cold protection insert, screwed tight ⚠ i.C. with wheel size 24 x 1 inch from RSH 390mm-430mm i.C. with wheel size 25 x 1 inch from RSH 420mm-460mm i.C. with wheel size 26 x 1 inch from RSH 450mm-510mm	1,6
<input type="checkbox"/> 100/2624	Side guard made of aluminum in frame color with cold protection insert, screwed tight ⚠ i.C. with wheel size 24 x 1 inch to HSH 380mm i.C. with wheel size 25 x 1 inch to HSH 410mm i.C. with wheel size 26 x 1 inch to HSH 440mm	1,6
<input type="checkbox"/> 783	Carbon arm support	1,6
<input type="checkbox"/> 748	Aluminum side guard in frame color with cold protection insert including height-adjustable armrest (screwed)	1,5,6
<input type="checkbox"/> 11	Side guard removable ⚠ only in combination with CODE 100, 783 oder 748	

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

13 Footrests

LSL **360** 370 380 390 400 410 420 430 440 450 460 470 480

CODE xR	LSL (mm)	DESCRIPTION
<input type="checkbox"/> 312	<input checked="" type="checkbox"/> A	Fußstütze durchgehend aus stabilem Aluminiumprofil in Rahmenfarbe hochklappbar, höhen und winkelverstellbar
<input type="checkbox"/> 311	<input checked="" type="checkbox"/> A	Footrest made of aluminum in frame color can be folded up, height and angle adjustable
<input type="checkbox"/> 4312	<input checked="" type="checkbox"/> A	LED footrest one-piece, made of strong aluminum in frame color, height / depth and angle adjustable, incl. Bluetooth module ⚠ i.C. with CODE 4959 from SW 360mm i.C. with CODE 4950 from SW 400mm
<input type="checkbox"/> 4313	<input checked="" type="checkbox"/> A	LED footrest one-piece, made of strong aluminum in frame color, height / depth and angle adjustable ⚠ i.C. with CODE 4959 from SW 360mm i.C. with CODE 4950 from SW 400mm
<input type="checkbox"/> 359		Calf belt

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

14 Frame colours

CODE	DESCRIPTION
<input type="checkbox"/> 130	Titanium grey
<input type="checkbox"/> 183	Metallic silver
<input type="checkbox"/> 186	Magic Blue
<input type="checkbox"/> 216	Deep black
<input type="checkbox"/> 206	Matt black
<input type="checkbox"/> 224	Blue-Metallic
<input type="checkbox"/> 227	Red Metallic
<input type="checkbox"/> 4023	Signal white
<input type="checkbox"/> 4040	Downtown Gray
<input type="checkbox"/> 4041	Magicblack Matt
<input type="checkbox"/> 4042	Blue Metallic Matt
<input type="checkbox"/> 4043	Brown metallic matt
<input type="checkbox"/> 4044	Red Metallic Matt
<input type="checkbox"/> 4045	Signal white matt
<input type="checkbox"/> 4046	Green Metallic Matt
<input type="checkbox"/> 4047	Purple metallic matt
<input type="checkbox"/> 4048	Copper metallic
<input type="checkbox"/> 4049	Olive green Matt
<input type="checkbox"/> 4913	Cross orange
<input type="checkbox"/> 4922	Metallic green




Other frame colors are available on request for an extra charge!

Please select by ticking the relevant boxes


*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

Pedal tubes

CODE	DESCRIPTION
<input type="checkbox"/> 309	Pedal caps (1 pair)  not i.C. with CODE 728/729/730 Support castors
<input type="checkbox"/> 4930	Pedal cap, mounted left  not i.C. with CODE 729/730 Support castors
<input type="checkbox"/> 4931	Pedal cap, mounted right  not i.C. with CODE 728/730 Support castors

Accessories

CODE	DESCRIPTION
<input type="checkbox"/> 833	Retaining strap with buckle
<input type="checkbox"/> 970	Walking aid holder  only in combination with CODE 309/4390/4391
<input type="checkbox"/> 819	Reflectors
<input type="checkbox"/> 951	Accessory net

Xtrends accessories

Item	DESCRIPTION
<input type="checkbox"/> 1080617	Tire pressure monitoring system for vehicles with 4 wheels
<input type="checkbox"/> 1080618	Tire pressure monitoring system for vehicles with 2 wheels

Remarks		
<div style="background-color: #e6f2ff; height: 250px;"></div>		
Date	City	Signature
<div style="background-color: #e6f2ff; height: 30px;"></div>	<div style="background-color: #e6f2ff; height: 30px;"></div>	<div style="background-color: #e6f2ff; height: 30px;"></div>