


Clou

9.500

 **Allowed user weight**
120 kg






Long-life LED lighting available as an option



Folded



30° back adjustment can be selected as option



Adjustable seat and back

The small folding one

EFFICIENCY

- Home and inpatient application
- Great life span through precision steel tube frame and high tensility of important components
- Optimal reinstallment through simple adaptation possibilities to the user
- Different seat widths in one frame
- Long term safety in spare part provision through a wide spread specialist dealer network
- Little need for repairs by implementing high value and durable parts

THERAPEUTIC BENEFIT

- Ideal for use in confined spaces due to compact dimensions
- Driving electronics adjustable in length and height, enabling comfortable, fatigue-free driving in an ergonomic posture
- Interactive programmable driving characteristics ensure individual adaptation
- Easy-to-operate control buttons enable handling even with very little strength
- Close driving to objects and easy transfer through leg supports removable to the side

FIELD OF USE

- For users that are no longer able to move around with a manual wheelchair, even indoors
- For extreme functional limitations of the lower extremities and little hand and arm strength
- For driving in confined spaces
- As transport aid for operation through an accompanying person

ORDER

COST ESTIMATE

CUSTOM DESIGN

Please fax the completed order form to: 05733 922 9311

Customer service: 05733 922 311

Billing address / Customer ref. no.:

Delivery address:

Consignment / Comments:

Order quantity:

MEYRA



Allowed user weight

120 kg



Hints

Our Terms and Conditions as stated at www.meyra.de/AGB shall apply
You can find information on the measurement process at www.meyra.de/messverfahren

Valid from 01/01/2021 - 12/31/2021.
Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check



Technical specifications

¹⁾ With standard legrest CODE 93/808

²⁾ under test conditions with drive battery 38 Ah (5h), 45 Ah (20h)

Seat widths in mm	380 - 500 / 430 - 550
Seat depth in mm	430 ¹⁾
Back support height in mm	400
Lower shank length in mm	360 - 500 ¹⁾
Arm support height in mm	190 - 300
Front seat height in mm	510
Length with leg supports in mm	1100 ¹⁾
Length without leg supports in mm	790
Width min./max. in mm	590 / 700
Width folded in mm	380
Height without head support in mm	910
Transport length in mm	790
Transport height min. in mm	910
Front wheel size in mm	200 x 50 (8")
Rear wheel size in mm	320 x 60 (12.5")
Max. obstacle height in mm	60
Turning radius in mm	1000

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

00 Chassis/Motorisation

CODE	SPEED	MOTOR POWER	RANGE
<input type="checkbox"/> 115	6 km/h	2 x 180 W	max. 30 km

01 Seat width (SW)

CODE	SW	DESCRIPTION
<input type="checkbox"/> 38	SW 380 - 500 mm	Seat width, continuously adjustable through the arm supports from 380 - 500 mm, (seat surface width 430 mm)
<input type="checkbox"/> 43	SW 430 - 550 mm	Seat width, continuously adjustable through the arm supports from 430 - 550 mm, (seat surface width 430 mm)

02 Seat system

CODE	DESCRIPTION
<input type="checkbox"/> 738	Standard seat belt without cushion

Seat system options

CODE	DESCRIPTION
<input type="checkbox"/> 560	Seat cushion 30 mm

03 Back systems

CODE	DESCRIPTION
<input type="checkbox"/> 737	Standard back, back support height 400 mm
<input type="checkbox"/> 736	Adjustable back, back support height 400 mm
<input type="checkbox"/> 26	30° Back adjustment through pneumatic spring

04 Arm supports and arm pads

CODE	DESCRIPTION
<input type="checkbox"/> 106	Standard arm support with upholstered arm pad and black clothes guard

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

05 Leg supports and footrests

CODE	DESCRIPTION
<input type="checkbox"/> 93/798	Legrest, lift off and swing away, divided footrest, hight adjustable, flip up, angle adjustable
<input type="checkbox"/> 92/798	Legrest, swing-up, lift off and swing away, divided footrests, hight adjustable, flip up and angle adjustable

06 Steering wheel tyres

CODE	DESCRIPTION
<input type="checkbox"/> 157/481	Pneumatic tyres 200 x 50 mm (8")
<input type="checkbox"/> 157/900	Puncture safe tyres 200 x 50 mm (8")

07 Tyres drive wheels

CODE	DESCRIPTION
<input type="checkbox"/> 815/482	Pneumatic tyres 20 x 60 mm (12.5")
<input type="checkbox"/> 815/904	Puncture safe tyres 20 x 60 mm (12.5")

08 Batteries and Chargers

CODE	DESCRIPTION
<input type="checkbox"/> 4696	Maintenance-free batteries 38Ah (5h) / 45 Ah (20h)
<input type="checkbox"/> 731	Battery charger 24 V / 6 A

09 Operating modules and brackets

CODE	DESCRIPTION
<input type="checkbox"/> 408	Operating module VR 2
<input type="checkbox"/> 60	Mounted left
<input type="checkbox"/> 74	Assembly right
<input type="checkbox"/> 4859	Operating module bracket, length adjustable
<input type="checkbox"/> 4860	Operating module bracket with parallelogram length and height adjustable

Please select by ticking the relevant boxes *SBI - Sonderbau/Individuallösung Standard version: If no other selection is made, the standard version will be delivered.

10 Lighting

CODE	DESCRIPTION
<input type="checkbox"/> 932	Active lighting

11 Frame- and cover colours

CODE	DESCRIPTION
<input type="checkbox"/> 4919	Silverline
<input type="checkbox"/> 287	Seat and back strap black

20 Accessories

CODE	DESCRIPTION
<input type="checkbox"/> 970	Walking aid holder
<input type="checkbox"/> 833	Retaining strap with buckle

Remarks		
<div style="background-color: #e0e0ff; height: 250px;"></div>		
Date	City	Signature
<div style="background-color: #e0e0ff; height: 30px;"></div>	<div style="background-color: #e0e0ff; height: 30px;"></div>	<div style="background-color: #e0e0ff; height: 30px;"></div>