



SMART S

2.370

 **Allowed user weight**
150 kg

 **Suitable for passenger transport**
Tested according to ISO 7176-19 (max. user weight 136 kg)





Optimised folding system
with compact folding
dimensions



Extremely low actuating forces



Individually adjustable
chassis



Unique side panel with
height adjustable armrest,
one handed operation

New Active Line

EFFICIENCY

- High part compatibility
- Little wear and proneness to repair through integral design of the frame and the quality of the attachment components
- Quick and long term assurance of spare parts provision through an area wide dealer structure

THERAPEUTIC BENEFIT

- Stabilization and mobilization of users with neurological diseases, through good support while sitting and low weight a.s.o.
- The straight front frame and swivel-away, removable leg supports make it easy to get close to everyday objects
- Folding mechanism for independent loading into a car
- Adaptable to functional changes of the user, e.g. by retrofitting electric drives
- Extensive equipment assortment for long-term users

FIELD OF USE

- For permanent and long term provision
- Very good suitability for many disease patterns, even for severely encumbered
- For users with remaining walking and standing abilities (standing transfer)
- With neurological diseases

ORDER

COST ESTIMATE

CUSTOM DESIGN

Please fax the completed order form to: 05733 922 9311

Customer service: 05733 922 311

Billing address / Customer ref. no.:

Delivery address:

Consignment / Comments:

Order quantity:

MEYRA



Allowed user weight

150 kg



Suitable for passenger transport

Tested according to ISO 7176-19 (max. user weight 136 kg)



Hints

Our Terms and Conditions as stated at www.meyra.de/AGB shall apply
You can find information on the measurement process at www.meyra.de/messverfahren

Valid from 01/01/2021 - 12/31/2021.
Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check



Technical specifications

¹⁾ When using auxiliary drives, the permissible total weight must not be exceeded

Seat widths in mm	320/340/360/380/400/420/440/460/480/500/520
Ready for use width in mm	Seat width + 180
Width folded in mm	280 with 0° wheel camber or 330 with 3° wheel camber
Seat depth in mm	380/400/420/440/460/480/500
Seat height front/rear in mm	400 to 540/370 to 530
Overall height in mm	Rear seat height + back support height
Back support height in mm	variable from 250 to 500 +25 mm
Arm support height in mm	ca. 170 - 245
Length with footrest in mm	880 - 1030 (depending on frame length, wheel size and axis position)
User weight	150 kg
Empty weight	from 10.5 kg
Permissible overall weight	165 kg ¹⁾
Transport weight	6,5 kg
Driving wheel	560mm (22") / 600 mm (24") / 635 mm (25") / 660 mm (26")
Brake test according to STANDARD	7°

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

00 Frame shape / seat width / seat depth

Seat width/Seat depth	Frame short 79°			Frame medium 79°		Frame large 79°	
	SD 380 mm	SD 400 mm	SD 420 mm	SD 440 mm	SD 460 mm	SD 480 mm	SD 500 mm
SW 320 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 340 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 360 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 380 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 400 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 420 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 440 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 460 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 480 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 500 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SW 520 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

including utensil pocket under the seat belt

Variable block

CODE	DESCRIPTION	TOTAL WIDTH
<input type="checkbox"/> 992	Varioblock for additional drives <i>! only necessary with selection of side panel CODE 81 or 84 in combination with Alber drives</i>	+ 20 mm

03 Castor wheel (without seat cushion) in combination with front seat height (FSH)

The maximum permissible seat height difference between the front and rear seat height is 90 mm

Code	Castor	Front seat height in mm (FSH)														
		400	410	420	430	440	450	460	470	480	490	500	510	520	530	540
151	Ø 100 Solid rubber	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
325	Ø 125 Solid rubber		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
327/908	Ø 142 Soft solid rubber black			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
327/909	Ø 142 Soft solid rubber grey			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
645	Ø 180 Pneumatic tires						<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
161	Ø 180 Solid rubber						<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Please select by ticking the relevant boxes

S *SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

X-Trend castor wheel (without seat cushion) in combination with front seat height (FSH)

Code	Castor	Front seat height in mm (FSH)											
		400	410	420	430	440	450	460	470	480	490	500	510
1068963	Ø 100 X-Trend Lightwheel	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1072385	Ø 100 X-Trend steering wheel	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1072386	Ø 125 X-Trend steering wheel		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Aluminum castor wheel (without seat cushion) in combination with front seat height (FSH)

FSH (mm)	400	410	420	430	440	450	460	470	480	490	500	510
A	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
B		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Code	Castor	Rim color					
		Silver	black	Green	Red	Blue	Gold
689	Ø 100 x 37 mm, Aluminum rim, Solid rubber	<input type="checkbox"/> A	<input type="checkbox"/> A	<input type="checkbox"/> A	<input type="checkbox"/> A	<input type="checkbox"/> A	<input type="checkbox"/> A
689	Ø 125 x 37 mm, Aluminum rim, Solid rubber	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B

03 Castors option

CODE	DESCRIPTION
<input type="checkbox"/> 903	Reinforced steering wheel forks, recommended when using additional drives and higher user weights

Castor connection

CODE	DESCRIPTION
<input type="checkbox"/> 2144	Reinforced steering head connection ! recommended when using additional drives

04 Seat cushions

CODE	DESCRIPTION
<input type="checkbox"/> 560	Seat cushion 30 mm
<input type="checkbox"/> 561	Seat cushion 60 mm
<input type="checkbox"/> 562	Velcro strip on seat belt

Please select by ticking the relevant boxes

S *SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

05 Rear seat height (RSH) in combination with drive wheel

The maximum permissible seat height difference between the front and rear seat height is 90 mm

RSH (mm)	370	380	390	400	410	420	430	440	450	460	470	480	490	500
A				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
B	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

CODE xR	RSH (mm)	DESCRIPTION	Ø
<input type="checkbox"/> 651	A	Standard wheel, incl. quick release axle (pair)	24"
<input type="checkbox"/> 2651	A	Standard wheel black, incl. quick release axle (pair)	24"
<input type="checkbox"/> 706	A	Drum braked wheel with radial spokes and drum brake for accompanying person ⚠ Only in conjunction with 0° wheel camber CODE 423 Only in conjunction with push handles CODE 412 or 502 Wheel covering CODE 98 is recommended.	24"
<input type="checkbox"/> 695	A	Active-wheel with radial spokes, high flange rim, hollow rim (Jumbo rim) in black ⚠	24"
<input type="checkbox"/> 752	A	Lightwheel incl. RightRun Light tyres and aluminium handrim black ⚠ only in combination with 0° wheel camber CODE 423 (pair) Recommended for user weight up to 90 kg	24"
<input type="checkbox"/> 848	A	SILVERDESIGN with radial spokes and high flange rim incl. quick-release axle (pair)	24"
<input type="checkbox"/> 651	B	Standard wheel, incl. quick release axle (pair)	22"
<input type="checkbox"/> 706	B	Drum braked wheel with radial spokes and drum brake for accompanying person ⚠ Only in conjunction with 0° wheel camber CODE 423 Only in conjunction with push handles CODE 412 or 502 Wheel covering CODE 98 is recommended.	22"

Please select by ticking the relevant boxes

S *SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

Drive wheels X-Trends

RSH (mm)	400	410	420	430	440	450	460	470	480	490	500	510	520	530
A	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
C		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
D			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
E				<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Item	RSH (mm)	DESCRIPTION	Ø
<input type="checkbox"/> 23700062	D	Day-Lite drive wheel black, including Speedlite tires and aluminum handrims black	25x1"
<input type="checkbox"/> 1068705	A	Six-Star incl. 8 bar high pressure tyres and aluminium handrim	24x1"
<input type="checkbox"/> 1073641	A	Spinergy-wheel 12 spokes black, incl. Marathon-Plus-Evolution and Tetra-handrim. Upon inquiry white or orange spokes possible!	24x1"
<input type="checkbox"/> 1073640	A	Spinergy-wheel 12 spokes black, incl. Marathon-Plus-Evolution and stainless steel handrim. Upon inquiry white or orange spokes possible	24x1"
<input type="checkbox"/> 1068713	A	Black Light incl. Marathon-Plus-Evolution and aluminium handrim	24x1"
<input type="checkbox"/> 1081354	A	Active wheel with radial spoke, Schwalbe ONE tires, hub, rim and handrims BLUE	24x1"
<input type="checkbox"/> 1081355	A	Active wheel with radial spoke, Schwalbe ONE tires, hub, rim and handrims GOLD	24x1"
<input type="checkbox"/> 1081356	A	Active wheel with radial spoke, Schwalbe ONE tires, hub, rim and handrims RED	24x1"
<input type="checkbox"/> 1068707	C	CROSS-wheel incl. aluminium rim and rough cleat profile tyres and stainless steel handrim	24x1.75"
<input type="checkbox"/> 1075478	D	Active-wheel with radial spokes, high flange hub, high flange hub, double walled rim (jumbo rim) in black, incl. Marathon-PLUS-Evolution and stainless steel handrim ⚠ Recommended from 100 kg user weight or in case of dynamic driving habits!	25x1"
<input type="checkbox"/> 1075479	E	Active-wheel with radial spokes, high flange hub, high flange hub, double walled rim (jumbo rim) in black, incl. Marathon-PLUS-Evolution and stainless steel handrim ⚠ Recommended from 100 kg user weight or in case of dynamic driving habits!	26x1"

⚠ not in conjunction with hand rims and tires

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

Horizontal drive wheel position (balance point setting)

	ACTIVE (Pos. - mm)			MEDIUM (Pos. - mm)		PASSIV (Pos. - mm)	
	1 - 100**	2 - 90**	3 - 80	4 - 70	5 - 60	6 - 50	7 - 40
Position	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

! Measurement: Centre of back tube to centre of axle in mm

Drive wheel position 1 and 2 cannot be fitted in conjunction with frame short CODE 351

Drive wheel options

CODE	DESCRIPTION
<input type="checkbox"/> 697	Release for tetra plegics, makes it easier to take off the drive wheel ! Only in conjunction with drive wheels CODE 848 oder CODE 695 recommended in conjunction with 3° wheel camber CODE 447

Wheel camber

CODE	DESCRIPTION
<input type="checkbox"/> 423	0° wheel camber (track adjustable, wheel camber can be altered subsequently)
<input type="checkbox"/> 447	3° wheel camber (track adjustable, wheel camber can be altered subsequent-ly) in combination with drum brake CODE 706 restricted seat height variation

Please select by ticking the relevant boxes

S *SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

07 Handrims and handrim options (not combinable with X-Trends wheels)

CODE	DESCRIPTION
<input type="checkbox"/> 445	Titanium, with 6 brackets ⚠ Only in combination with active wheel CODE 695,
<input type="checkbox"/> 493	Stainless steel, polished
<input type="checkbox"/> 494	Aluminium ⚠ Not in combination with active wheel CODE 695
<input type="checkbox"/> 168	Optigripp-special handrim for tetraplegics
<input type="checkbox"/> 849	Curve L, aluminium black anodised, ergonomic handrim profile, H 31 mm ⚠ only in combination wheel CODE 695, 651 or 2651
<input type="checkbox"/> 850	Curve Grip L, anti-slip-coated, aluminium, suited for tetraplegics ergonomic handrim profile, H 31 mm ⚠ only in combination wheel CODE 695, 651 or 2651
<input type="checkbox"/> 853	Curve GEKKO, aluminum black anodized with exchangeable rubber lip, ergonomic handrim profile, H 31 mm ⚠ only in combination wheel CODE 695, 651 or 2651
<input type="checkbox"/> 871	Double handrim for seat widths up to 480 mm ⚠ Only in conjunction with PU tyres CODE 902 Not in conjunction with drum brakes CODE 706 not suitable for very active use
<input type="checkbox"/> 166	Handrim coating made of silicone, smooth, black ⚠ Not in combination with active wheel CODE 695
<input type="checkbox"/> 692	Handrim distance narrow
<input type="checkbox"/> 693	Handrim distance far

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

08 Tyres (not combinable with X-Trends wheels)

CODE	DESCRIPTION	Ø
<input type="checkbox"/> 486	Pneumatic tires 4 bar	22x1 3/8" - 24x1 3/8"
<input type="checkbox"/> 496	Pneumatic tires 8 bar	22x1" - 24x1"
<input type="checkbox"/> 902	PU-tyres, puncture safe ⚠ Not in combination with active wheel CODE 695	22x1" - 24x1"
<input type="checkbox"/> 902	PU-tyres, puncture safe ⚠ Not in combination with active wheel CODE 695	22x1 3/8" - 24x1 3/8"
<input type="checkbox"/> 897	Marathon-Plus-Evolution tyres, puncture protected	22x1" - 24x1"
<input type="checkbox"/> 363	Schwalbe ONE tires, no profile, puncture resistant	22x1" - 24x1"
<input type="checkbox"/> 395	Maxxis VN Master Alpha Protect tire black, puncture resistant, 10 bar	22x1" - 24x1"

⚠ Tyres consist of a rubber compound which may leave permanent marks or be hard to remove from some surfaces (e.g. plastic, wooden or parquet flooring, rugs and carpets).
 ☒ We cannot accept liability for any abrasion damage to flooring.

09 Brakes

CODE	DESCRIPTION
<input type="checkbox"/> 305	Pressure brake for user
<input type="checkbox"/> 60/650	One-hand pressure brake left (incl. Brake lever long, hinged)
<input type="checkbox"/> 74/650	One-hand pressure brake right (including brake lever long, hinged)
<input type="checkbox"/> 73	Side panel brake, brake integrated in clothing guard ⚠ only in conjunction with CODE 100
<input type="checkbox"/> 681	Brake lever extension, foldable, firmly mounted

10 Back strap height (BH)

DESCRIPTION	Back strap height in mm					
	250	300	350	400	450	500
Back strap, adaptable via Velcro straps adjustable +25 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Please select by ticking the relevant boxes

S *SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

11 Back support

CODE	DESCRIPTION
<input type="checkbox"/> 410	Back with push handle, short
<input type="checkbox"/> 412	Back with push handle, long
<input type="checkbox"/> 141	Push handles, swivel-away
<input type="checkbox"/> 502	Back support with push handle mounted to the frame tube, continuously height adjustable to 150 mm ⚠ Not in combination with push handle code 140/141
<input type="checkbox"/> 140	Push handle integrated in the back tube, height adjustable up to 100 mm ⚠ Not in combination with lumbar bent code 622

Back support options









CODE	DESCRIPTION
<input type="checkbox"/> 29	Back, angle-adjustable from seat surface, +/-15° in 5° increments, can be adjusted subsequently ⚠ from back height 300 mm Only in conjunction with push handle CODE 410/ 412. Not in conjunction with lumbar contour CODE 622 not i.c.w. back height 250 mm
<input type="checkbox"/> 621	Back without lumbar support
<input type="checkbox"/> 622	Back with lumbar support ⚠ possible from BH 300 mm
<input type="checkbox"/> 2032	ErgoSeat Pro thorax pelottes, small, pair ⚠ Recommended at back height 300 - 400 mm
<input type="checkbox"/> 2021	ErgoSeat Pro thorax pelottes, large, pair ⚠ Recommended at back height 400 - 450 mm
<input type="checkbox"/> 814	Horizontal push bar, stabilises the back support under heavy strain ⚠ not in combination with push handles CODE 141. Recommended from SB 460
<input type="checkbox"/> 913	Head support, 300 mm continuously height adjustable ⚠ only in combination with push bar CODE 814

Please select by ticking the relevant boxes


*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

12 Arm supports and arm pads

CODE	DESCRIPTION	Functions
<input type="checkbox"/> 100	Aluminium arm support in frame colour with carbon insert  from rear seat height 440 mm	1,6
<input type="checkbox"/> 624/100	Aluminium arm support in frame colour with carbon insert  to rear seat height 430 mm	1,6
<input type="checkbox"/> 98	Side part made of aluminum in frame color with carbon insert and 50 mm wheel cover  from rear seat height 440 mm	1,6
<input type="checkbox"/> 624/98	Side part made of aluminum in frame color with carbon insert and 50 mm wheel cover  to rear seat height 430 mm	1,6
<input type="checkbox"/> 81	Side guard, armrest short, height adjustable with one hand operation  from rear seat height 440 mm	2,3,4,5
<input type="checkbox"/> 84	Side guard, armrest long, height adjustable with one hand operation  from rear seat height 440 mm	2,3,4,5
<input type="checkbox"/> 748	Aluminium arm support in frame colour with carbon insert, incl. height adjustable arm support  Not in combination with light brake Code 595 from rear seat height 440 mm	1,5,6
<input type="checkbox"/> 783	Carbon arm support	1,6
<input type="checkbox"/> 302	Aluminium Light arm support in frame colour  from the rear seat height 430 mm	1

12 Side panel option


CODE	DESCRIPTION
<input type="checkbox"/> 24	Swing-away side panel bracket  only in combination with Code 98, 100 783 or without side guard



Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

13 Leg supports and footrests

CODE	DESCRIPTION
92	Leg supports, swivelling upward, removable and swivelling
 from seat width 380mm	
93	Legrest upper part standard, detachable and swing-away, front frame angle 70°
293/3532	like Code 93 with shortened lower leg length
293	Legrest upper part compact, detachable and swing-away, front frame angle 78°
94	Without leg supports and footrest
4507/759	amputee leg support, single, right-mounted only
4507/760	Amputee leg support, single, left-mounted only

798	Plastic footrest, split, fold-up, height-adjustable and angle-adjustable
808	Single foot plates made of aluminum in frame color, can be folded up, with depth and angle adjustment
312	Footrest made of sturdy aluminum profile in frame color, foldable, height and angle adjustable
315	Footrest made of aluminum in frame color can be folded up, height and angle adjustable
4312	Stable aluminium profile LED footrest incl. headlamp and tail light in frame colour, one-piece, height, depth and angle-adjustable, incl. Bluetooth module (from seat width 360 mm up)
 From seat width 360 mm	
4313	Stable aluminium profile LED footrest incl. headlamp and tail light in frame colour, one-piece, height, depth and angle-adjustable (from SW 360 mm up)
 From seat width 360 mm	

Please select by ticking the relevant boxes

S *SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

Leg supports and footrests

The meaning of the individual codings can be found in the previous section

LSL (mm)	270	280	290	300	310	320	330	340	350	360	390	400	410	420	430	440	450	460	470	480	490	
A	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>												
B											<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
C												<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

leg support	Footrests CODE (USL)						
	798	808	312	315	4312	4313	94
92	<input type="checkbox"/> C						
93		<input type="checkbox"/> C	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B	
93/3532		<input type="checkbox"/> A	<input type="checkbox"/> A	<input type="checkbox"/> A	<input type="checkbox"/> A	<input type="checkbox"/> A	
293		<input type="checkbox"/> C	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B	<input type="checkbox"/> B	
4507/759	<input type="checkbox"/> C	<input type="checkbox"/> C					
4507/760	<input type="checkbox"/> C	<input type="checkbox"/> C					
94							<input type="checkbox"/>

! CODE 4312 only in connection from seat width 360 mm up
 CODE 4313 only in connection from seat width 360 mm up

Footrest accessories

CODE	DESCRIPTION
<input type="checkbox"/> 359	Calf belt
<input type="checkbox"/> 925	Operating aid for locking device, in case of limited hand function ! only in combination with upper leg support compact CODE 293

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

14 Frame colours

CODE	DESCRIPTION
<input type="checkbox"/> 130	Titanium grey
<input type="checkbox"/> 227	Redmetallic
<input type="checkbox"/> 183	Metallic silver
<input type="checkbox"/> 188	Blackstar
<input type="checkbox"/> 216	Deep black
<input type="checkbox"/> 4023	Signal white
<input type="checkbox"/> 186	Magic Blue
<input type="checkbox"/> 224	Blue-Metallic
<input type="checkbox"/> 206	Matt black
<input type="checkbox"/> 4040	Downtowngrey
<input type="checkbox"/> 4041	Magicblack Matt
<input type="checkbox"/> 4042	Blue Metallic Matt
<input type="checkbox"/> 4043	Brown metallic matt
<input type="checkbox"/> 4044	Red Metallic Matt
<input type="checkbox"/> 4045	Signal white matt
<input type="checkbox"/> 4046	Green Metallic Matt
<input type="checkbox"/> 4047	Purple metallic matt

15 Back strap colours (seat strap always in black)

CODE	DESCRIPTION
<input type="checkbox"/> 235	Trevira black
<input type="checkbox"/> 292	Net fabric, breathable , black
<input type="checkbox"/> 296	Net fabric, breathable , red

16 Wheelchair accessories

Hand protectors and spoke guards

CODE	DESCRIPTION
<input type="checkbox"/> 489	Transparent

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

Support castors

CODE	DESCRIPTION
<input type="checkbox"/> 691	Support castors, 1 pair, swivelling (recommended) ⚠ not in combination with step tube CODE 309
<input type="checkbox"/> 712	1 support castor, swivelling, mounted left ⚠ Not in combination with pedal tube CODE 309/4930, not in combination with shunting castors CODE 746
<input type="checkbox"/> 713	1 support castor, swivelling, mounted right ⚠ Not in combination with pedal tube CODE 309/4931, not in combination with shunting castors CODE 746

Accessories

CODE	DESCRIPTION
<input type="checkbox"/> 746	1 pair shunting castors, e.g. for tight spaces or in an aeroplane ⚠ Not i.c.w. pedal tube CODE 309 or support castors CODE 691, from RSH 420 mm
<input type="checkbox"/> 833	Retaining strap with buckle
<input type="checkbox"/> 970	Walking aid holder ⚠ *nur i.V.mit Auftrittrohr CODE 309/4930/4931 oder CODE 691/712/713
<input type="checkbox"/> 819	Reflectors
<input type="checkbox"/> 951	Accessory net
<input type="checkbox"/> 309	Step tubes (1 pair)
<input type="checkbox"/> 4930	1 step tube, mounted on the left ⚠ not in combination with anti-tip-wheels, CODE 691/712 and anti-tip-wheels, CODE 746
<input type="checkbox"/> 4931	1 step tube, mounted on the right ⚠ not in combination with anti-tip-wheels, CODE 691/713 and anti-tip-wheels, CODE 746
<input type="checkbox"/>	205542100 Hand compressor
<input type="checkbox"/> 679	Folding lock
<input type="checkbox"/> 929	Therapy table plastic, transparent, adjustable from 380 - 480 mm seat width ⚠ only in combination with side panel CODE 84
<input type="checkbox"/> 930	Therapeutic table for manual wheelchairs made of plastic, transparent, for sliding on, for defined seat width ⚠ only in combination with side panel CODE 84

⚠ With 24" wheels and a rear seat height of less than 470 mm, there is a resulting limitation in terms of armrest height with use of the therapy table Code 929/930

Please select by ticking the relevant boxes *SBI - Sonderbau/Individuallösung Standard version: If no other selection is made, the standard version will be delivered.**Xtrends accessories**

Item	DESCRIPTION	HMV
<input type="checkbox"/> 261519300	X-TREME Backpack	
<input type="checkbox"/> 261519600	Wheel pouch „BLACK WHEEL CASE“ fitting for max. 2 x 24“ wheels resp. each 1 x 25“ or 26“ wheel	18.99.99.0999
<input type="checkbox"/> 1074231	Wheel pouch „MEYRA ACTIVE“, fitting for max. 3 x 28“ wheels	
<input type="checkbox"/> 206964600	3M Spoke reflectors	18.99.99.0999
<input type="checkbox"/> 1080617	Tire pressure monitoring system for vehicles with 4 wheels	
<input type="checkbox"/> 1080618	Tire pressure monitoring system for vehicles with 2 wheels	

Remarks		
<div style="background-color: #e6f2ff; height: 250px;"></div>		
Date	City	Signature
<div style="background-color: #e6f2ff; height: 30px;"></div>	<div style="background-color: #e6f2ff; height: 30px;"></div>	<div style="background-color: #e6f2ff; height: 30px;"></div>