


Budget

9.050 Stock version

 **Allowed user weight**
130 kg

 **Suitable for passenger transport**
Tested to ISO 7176-19





Easy removal function of the side plate



Knee lever brake



Lockable armrests can be folded back



Front seat height adjustable without changing parts

Economic and innovative

EFFICIENCY

- Specially designed for the current provision situation
- Technical variability through compatibility with EUROCHAIR: Arm supports, leg supports, seat- and back straps, pressure and drum brake

THERAPEUTIC BENEFIT

- Easy to fold for easy handling in everyday life

FIELD OF USE

- Push- handrim drive
- Home and inpatient application
- Rheumatism
- Geriatric provision

ORDER

COST ESTIMATE

CUSTOM DESIGN

Please fax the completed order form to: 05733 922 9311

Customer service: 05733 922 311

Billing address / Customer ref. no.:

Delivery address:

Consignment / Comments:



Basic equipment

- Drive wheel 24" aluminium
- Steering wheel 180 x 50 mm
- Puncture-proof tyres
- Arm supports up-foldable
- Leg supports, removable and swivelling
- Seat height 510 mm

- Seat depth 430 mm
- Back height 420 mm
- Handrim distance narrow
- Frame colour Silverline
- Back with lumbar support



Allowed user weight

130 kg



Suitable for passenger transport

Tested to ISO 7176-19



Hints

Our Terms and Conditions as stated at www.meyra.de/AGB shall apply

When transported in public transport (e.g. train), the dimensions of the wheelchair may be higher than the guide values of the transport company. Please inquire in advance.

*Bei Verwendung von Zusatzantrieben darf das zul. Gesamtgewicht nicht überschritten werden.

Valid from 01/01/2021 - 12/31/2021.

Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de

You avoid queries and incorrect orders thanks to the buildability check



Technical specifications

Seat widths in mm	380 / 400 / 430 / 460 / 480 / 510
Ready for use width in mm	seat width + 180
Width folded in mm	280
Seat depth in mm	430
Front seat height in mm	470 - 510
Seat height in the back in mm	440 - 490
Overall height in mm	960
Back support height in mm	420
Arm support height in mm	210,5
Length with footrest in mm	1040
User weight	130 kg
Empty weight	18,5 kg
Permissible overall weight	148,5 kg
Driving wheel	610 mm (24")

Please enter your quantities in the boxes

01 Pressure Brake

Item	DESCRIPTION
<input type="text"/> 9905141	SW 380 Pressure Brake
<input type="text"/> 9905142	SW 400 Pressure Brake
<input type="text"/> 9905143	SW 430 Pressure Brake
<input type="text"/> 9905144	SW 460 Pressure Brake
<input type="text"/> 9905145	SW 480 mm Pressure Brake
<input type="text"/> 9905146	SW 510 mm Pressure Brake

02 Drum Brake

Item	DESCRIPTION
<input type="text"/> 9905147	SW 380 Drum Brake
<input type="text"/> 9905148	SW 400 Drum Brake
<input type="text"/> 9905149	SW 430 Drum Brake
<input type="text"/> 9905150	SW 460 Drum Brake
<input type="text"/> 9905151	SW 480 Drum Brake
<input type="text"/> 9905152	SW 510 Drum Brake

Please enter your quantities in the boxes

20 Accessory (is sold separately)

Item	DESCRIPTION
<input type="text"/>	10622120092 Upper leg support Code 92 left
<input type="text"/>	10622130092 Upper leg support Code 92 right
<input type="text"/>	10575320092 Upper leg supports Code 92 (pair)
<input type="text"/>	1045844 Brake lever extension Code 986
<input type="text"/>	1063168 Walking aid holder
<input type="text"/>	205651500 Half-sided therapy table left Code 931
<input type="text"/>	205651600 Half-sided therapy table right Code 931
<input type="text"/>	1041414 Retaining strap with buckle
<input type="text"/>	1072681 Height adjustable push handles Code 502 (refit kit)
<input type="text"/>	8110833 Knee cushion black for Code 92 (pair)
<input type="text"/>	1032429 Air-pump with holder
<input type="text"/>	1047486 Reflectors
<input type="text"/>	1041414 Safety belt code 833
<input type="text"/>	1069874 Safety package, consisting of: Retaining strap, cane holder, passive lighting, support castors
<input type="text"/>	10598250691 Support castors, length adjustable, removable (1 pair)
<input type="text"/>	1064946 Therapy table plastic, transparent, adjustable from 38 - 48 cm seat width
<input type="text"/>	1064946 Therapy table Code 929 SW 38 to 48 cm adjustable
<input type="text"/>	1064947 Therapy table Code 929 SW 48 to 58 cm adjustable
<input type="text"/>	205642506 Universal backpack Code 962
<input type="text"/>	205538400 ACCESSORY NET CODE 951

Please enter your quantities in the boxes

30 **Seat cushion (supplied separately)**

	Item	SW-SD-CH (mm)	DESCRIPTION
<input type="text"/>	1006685	380-430-30	Seat cushion black
<input type="text"/>	1006685	400-430-30	Seat cushion black
<input type="text"/>	1006825	430-430-30	Seat cushion black
<input type="text"/>	1015928	460-430-30	Seat cushion black
<input type="text"/>	1046382	480-430-30	Seat cushion black
<input type="text"/>	1046382	510-430-30	Seat cushion black

Remarks

Date

City

Signature