

Solero Light

9.072 Individual



Allowed user weight

130 kg up to seat width 460 mm

150 kg from seat width 480 mm



Suitable for passenger transport

Tested according to ISO 7176-19 (max. user weight 136 kg)





Adaptable in every respect



Simplicity in handling

EFFICIENCY

- Outstanding price-performance ratio
- Economic advantages through the compatibility to EUROCHAIR
- Puncture-proof tyres

THERAPEUTIC BENEFIT

- Anatomically formed seat- and back cushions ensure sitting comfort
- Seat- and back settings can be continuously adjusted
- Head support is height-, depth- and angle adjustable
- Very slim overall width has advantages when driving through narrow doors or in packed rooms
- Ergonomically shaped push bar. Facilitates pushing in tilted position and allows different grip options
- Pedal tube and special gas pressure suspension provide an especially gentle adjustment with a minimum of strength required.

FIELD OF USE

- Geriatric care in senior and care centers as well as in the home
- For physical limitations, that require a supported sitting position
- Good suitability and mobilisation for self-drivers
- Userweight 130 to 150 kg

ORDER

COST ESTIMATE

CUSTOM DESIGN

Please fax the completed order form to: **05733 922 9311**

Customer service: **05733 922 311**

Billing address / Customer ref. no.:

Delivery address:

Consignment / Comments:

Order quantity:

ORTOPEDIA
BY MEYRA



Basic equipment

- Seat angle from -2 up to +25°
- Back angle adjustable from 0 up to 45°
- Ergonomical Seat- and Back pads
- Adjustable seat-depth from 420 up to 480 mm
- Headrest
- Puncture-proof tyres
- Drive wheels 24"
- Steering wheel 200 x 50 mm



Allowed user weight

130 kg up to seat width 460 mm
150 kg from seat width 480 mm



Suitable for passenger transport

Tested according to ISO 7176-19 (max. user weight 136 kg)

- Handrim aluminium
- Detachable Side part
- Armrest height- and depth-adjustable
- leg support detachable, swivelling away-, and flip-up
- Footrests flip-up, height and angle-adjustable
- Frame colour Silverline
- Incontinence material black
- anti-tip wheels

**Technical specifications**

Seat widths in mm	430 - 460 / 480 - 510
Ready for use width in mm	min. 620
Seat depth in mm	420 - 480
Seat height with seat padding in mm	440 (520)
Overall height in mm	1040 - 1180
Back support height in mm	530 - 560
Back angle	0° - 45°
Seat camber	-2° - 25°
Arm support height in mm	210 - 270
Length with footrest in mm	1180
Length without footrests in mm	910
User weight	130 kg / 150 kg
Empty weight	38 kg
Permissible overall weight	168 kg / 198 kg
Transport weight	23 kg
Driving wheel	610 mm(24 inch)

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

01 Seat system, seat width (SW) and seat depth (SD)

CODE	DESCRIPTION	USER WEIGHT	SEAT DEPTH
<input type="checkbox"/> 43	SW 430 mm	max. 130 kg	Seat depth - modifiable from 420 to 480mm
<input type="checkbox"/> 46	SW 460 mm	max. 130 kg	Seat depth - modifiable from 420 to 480mm
<input type="checkbox"/> 48	SW 480 mm	max. 160 kg	Seat depth - modifiable from 420 to 480mm
<input type="checkbox"/> 51	SW 510 mm	max. 160 kg	Seat depth - modifiable from 420 to 480mm

Upholstered seat

CODE	DESCRIPTION
<input type="checkbox"/> 562	Removable, ergonomically formed seat cushion (washable cover) with seat plate

Seat system options

CODE	DESCRIPTION
<input type="checkbox"/> 4724	Abduction wedge
<input type="checkbox"/> 959	Cushion to reduce seat width by 50 mm
<input type="checkbox"/> 4613	Retaining strap with buckle

02 Back system and Back support height (BH)

CODE	DESCRIPTION	BH
<input type="checkbox"/> 110	Back padding anatomically shaped, Back angle adjustable from 0 - 45°	540 mm, can be altered to 570 mm

Back support options

CODE	DESCRIPTION
<input type="checkbox"/> 814	Horizontal push bar ergonomically formed

Back support accessories

CODE	DESCRIPTION
<input type="checkbox"/> 4739	Torso cushion, can be mounted both sides

Please select by ticking the relevant boxes

*SRL - Sonderbau/Individuallösung

04 Steering- /drive wheel in combination with brake

CODE	Driving wheel	Castor
		PU tyres puncture-proof 200 x 50 mm
815/641	12 x 1 1/12" firmly screwed on, 2 x pressure brake user	<input type="checkbox"/>
815/706	12 x 1 1/12" quick release axle, drum brake accompanying person, 2 x pressure brake user	<input type="checkbox"/>
896/651	24 x 1 3/8" 2 x pressure brake user	<input type="checkbox"/>
896/706	24 x 1 3/8" drum brake accompanying person, 2 x pressure brake user	<input type="checkbox"/>

05 Tyres

CODE	DESCRIPTION
<input type="checkbox"/> 902	PU-tyres, puncture safe

06 Handrims

CODE	DESCRIPTION
<input type="checkbox"/> 494	Aluminium

Handrims distance

CODE	DESCRIPTION
<input type="checkbox"/> 692	Handrim distance narrow (extendable to far)

07 Arm supports

CODE	DESCRIPTION	Functions
<input type="checkbox"/> 106	Upholstered arm support, removable, depth- and height adjustable	1,2,3
<input type="checkbox"/> 4426	Arm support transparent with height adjustable padded arm support	1,2,3

Please select by ticking the relevant boxes

*SBI - Sonderbau/Individuallösung

Standard version: If no other selection is made, the standard version will be delivered.

08 Leg supports

CODE DESCRIPTION

92 Leg supports, swivelling upward, removable and swivelling

09 Footrests

LSL (mm) **350 360 370 380 390 400 410 420 430 440 450 460 470 480 490 500 510**

CODE xR LSL (mm) DESCRIPTION depth

798 Footrests divided, height adjustable, up-folding and angle adjustable 160 mm

2092 steples adjustment of the lower leg length

10 Frame colours

CODE DESCRIPTION

4919 Silverline

11 Incontinency cover/colour

CODE DESCRIPTION

4980 Incontinency cover, black

12 Hand and spoke guard

CODE DESCRIPTION

4610 Transparent for 24"



13 Support castors

CODE DESCRIPTION

691 Safety support (pair)

Please select by ticking the relevant boxes *SBI - Sonderbau/Individuallösung Standard version: If no other selection is made, the standard version will be delivered.

14 Therapy tables

CODE	DESCRIPTION
<input type="checkbox"/> 929	Therapy table plastic, transparent, adjustable from 380 - 580 mm seat width  not i.c. with seat width 510mm CODE 51 and padded arm support CODE 106
 With 24" wheels and a rear seat height of less than 470 mm, there is a resulting limitation in terms of armrest height with use of the therapy table as the lowest armrest position is then not possible.	

15 Accessories

CODE	DESCRIPTION
<input type="checkbox"/> 819	Passive lighting (1 set)
<input type="checkbox"/> 4616	Walking aid holder
<input type="checkbox"/> 585	Infusion post

Remarks		
<div style="background-color: #e6f2ff; height: 250px;"></div>		
Date	City	Signature
<div style="background-color: #e6f2ff; height: 30px;"></div>	<div style="background-color: #e6f2ff; height: 30px;"></div>	<div style="background-color: #e6f2ff; height: 30px;"></div>